

# Airstream E Class II Microbiological Safety Cabinet; 1.2m

When safety is paramount!  
Operator protection, product protection and environment protection!



Package Offer

**Only \$10,970**

Usually \$15,673

Includes:

UV light and 2 power sockets

Code: ESCOAC2-4E8

*Pricing valid for current stock only*

**In Stock**

Optional adjustable height stand  
available for  
only **\$895** (usually \$1,515)  
just add code ESCOSTC-4AO



Sentinel Gold microprocessor controller  
with digital display shows cabinet status

- UV resistant glass side panels allow entry of “natural” light
- Manual sash automatically activates internal lighting and airflow when set to the marked working height
- Night set back mode reduces power consumption by a further 60%
- Ergonomic sloping front combined with tempered glass side panels allows entry of more natural light
- ULPA downflow and exhaust filters with 99.999% efficiency, 10X better than HEPA filter (99.99%) improves safety and delivers ISO Class 3 air cleanliness
- Raised armrest external to cabinet for comfortable working position
- Multi-piece stainless steel worksurface is easy to remove for cleaning
- Certified by PHE (Public Health England) for conformance to EN 12469
- Airflow sensor monitors real-time airflow for safety and alerts the user if airflow is insufficient
- Removable easy-to-clean paper catch
- Energy efficient DC ECM Motor ( with 70% energy savings compared to AC motor)
- ISOCIDE™ Powder Coat – painted surfaces have silver ion impregnated coat to inhibit microbial growth
- 3-year warranty