

# SENSIS*Strip* Lateral Flow Tests



**Only \$388**

Usually \$420


(except Gluten \$408)



see page 2 for list of tests available

The lateral flow assays are the perfect tool to control the potential cross-contamination of food allergens, avoiding external analyses and thus considerably reducing the time needed to release a food product batch to the market.

- High sensitivity, low limits of detection
- Semi-quantitative result provided by an evaluation card
- No need for a highly qualified end user: fast and easy to use
- All necessary & ready-to-use materials provided



**For qualitative, consistent & auditable results:** the RapidScan ST5-W lateral flow reader is unique to the market with its flexibility, multiple target analysis and its compatibility with applications for our allergen, GMO, glyphosate and algal toxins lateral flow tests.

This portable instrument performs well in both laboratories and in field testing environments.



Contact us for more details



Ref. Code:  
GOLD-SENSI

VALID UNTIL:  
31 July 2024



0800 34 24 66

sales.nz@bio-strategy.com

www.bio-strategy.com

shop.bio-strategy.com

# SENSIStrip Lateral Flow Tests

**Contains:** 20 pcs, Food matrix sensitivity: 2.5 ppm, Swabbing sensitivity: 0.007 µg/cm<sup>2</sup>,  
Rinse water sensitivity: 0.17 mg/L"

Code	Test
GOLDHU0030110	Almond
GOLDHU0030112	Casein
GOLDHU0030113	Cashew
GOLDHU0030114	Coconut
GOLDHU0030115	Crustaceans (Tropomyosin)
GOLDHU0030116	Egg White
GOLDHU0030119	Hazelnut
GOLDHU0030122	Macadamia Nut
GOLDHU0030123	Milk Protein
GOLDHU0030125	Mustard
GOLDHU0030128	Peanut
GOLDHU0030129	Pecan
GOLDHU0030130	Pistachio
GOLDHU0030131	Sesame
GOLDHU0030132	Total Soy
GOLDHU0030133	Walnut
GOLDHU0030118	Gluten* * <b>Only \$408</b> Usually \$440
<b>NEW</b>	
GOLDHU0030111	Beta-Lactoglobulin
GOLDHU0030117	Fish
GOLDHU0030120	Lupin
GOLDHU0030121	Lysozyme
GOLDHU0030126	Ovalbumin
GOLDHU0030127	Brazil Nut



Ref. Code:  
GOLD-SENSI

VALID UNTIL:  
31 July 2024



0800 34 24 66

sales.nz@bio-strategy.com

www.bio-strategy.com

shop.bio-strategy.com