



# A bestseller in the laboratory!

## RW 20 Digital Overhead Stirrer Package



Only \$7,995

Usually \$9,765

Code: IKAA9011605

### Package includes

- RW 20 digital
- Plate stand
- Boss head clamp
- Propeller stirrer, 4-bladed

- Robust, slimline, ergonomic design with digital display
- Constant power-drive, with two speed ranges from 60 – 2,000 rpm
- Push-through agitator shafts (only when stationary)
- For stirring quantities of up to 20L (H2O)



Wide range of other stirring elements and stands also available.  
Contact us for pricing.



Part of DKSH Group

Ref. Code:  
IKA-RW20

VALID UNTIL:  
31 July 2024



0800 34 24 66

sales.nz@bio-strategy.com

www.bio-strategy.com

shop.bio-strategy.com