


Aquastab Algaecide For Water Baths

Reduce water bath contamination with
Aquastab Algaecide!



Only \$75

Usually \$103

Code:SCAR0000LZB929

100ml Bottle

Aquastab Algaecide is a highly effective blue additive to reduce the contamination from bacteria and fungi in your water bath.

- Recommended temperature range 5°C-90°C
- Highly effective algaecide only to be used with water in bath thermostats.
- Composition: 2% Poly(dimethyl imino)-2-hydroxy-1,3-propanediyl chloride
- Use 5ml/10L of water - 100ml bottle makes up 200L

Marketed by Scharlau

[Download](#) Safety Data Sheet