

Isotherm[®] Forced Convection Laboratory Oven

Ambient + 7.5°C to 250°C
Supplied with 2 shelves

In Stock



watch >

170 Litre

Only \$4,775

Usually \$6,613

Code: ESCOOPA-170-8

240 Litre

Only \$5,975


Usually \$8,272

Code: ESCOOPA-240-8


These forced convection ovens are used for high-volume thermal convection applications.

They provide uniform temperature throughout the chamber necessary for annealing, drying, sterilizing, and other lab functions.

- Dual display showing set and actual temperature
- Isocide™ Antimicrobial Powder Coating: Eliminates 99.9% of surface bacteria within 24 hours of exposure
- Forced air circulation with adjustable fan speed
- Permanently lubricated maintenance-free fan
- Air flow slider to adjust exchange rate of air
- Microprocessor PID control with 10 program memory



Single-piece stainless steel chamber with directly mounted shelves making oven easy to clean



4-zone heated air jacket for maximum temperature uniformity

Prices and discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices in local currency excluding GST and local delivery charges. Products advertised may be for research use only (RUO) or in vitro diagnostic use (IVD). Please refer to manufacturer's product insert for further information. For product use outside of that specified on the manufacturer's labelling and/or instructions for use, Bio-Strategy does not take responsibility.

Ref. Code:
ESCO-OFA

VALID UNTIL:
31 July 2024