

Now \$16,785

Usually \$24,245

Code: MDF-DU502VX-PE



The secure ultra-low temperature freezers for the storage of high value samples

Offering 528 litres of storage safety due to the dual cooling system.

Two independent refrigeration systems provide the best redundancy you can have! – if one refrigeration circuit fails, the other circuit will maintain -70°C.



VIP Plus insulation: provides typically 24% more storage capacity compared to conventionally insulated freezers with the same footprint.

- Quiet in operation – only 52dB(A)
- Filter-less condenser design: eliminates the need for routine filter cleaning
- Automatic vacuum relief port for quicker access to samples after door closure
- Hot gas bypass combined with fivefold door gasket assists in minimising frosting and ice build up around door seal
- Three access ports
- Intelligent ECO Mode: optimising compressor operation to match the freezer's running condition and safe electricity costs
- User friendly colour graphical control panel with digital display
- Alarms for high & low temperature, door open and power failure
- USB port for downloading logged data
- Castors for mobility & locating feet when in working position
- Keyed door locks for security
- Dry remote alarm contacts
- Suitable for operation on 220/240V 50Hz 10amp supply
- 5-year warranty