

Stable Temperature



**Only 2
Available**

Now \$2,240

Usually \$3,210

Code: SANYLPR-400-PE

Includes 5 adjustable shelves



- Temperature Range: 4°C to 14°C selectable
- Cyclic Defrost (timer controlled)
- Forced air circulation ensures improved temperature uniformity
- Visual & audible alarms notify end user of temperature deviation outside "safe" range
- Two layer glass window in door for easy sample viewing
- Low noise design; 48dB(A)
- Access port on rear wall of cabinet for remote probe/cable access
- Outer door has cylinder keylock for security
- Suitable for operation on 240V 50Hz 10 amp supply
- 2-year Warranty



Easy to read display with digital microprocessor control