



# Disperser T 25

## Digital Ultra-Turrax Package

**Just \$6,750**

**Usually \$8,170**

Code: IKAA9024205

- Motor protection against overload
- Quick-connect coupling to exchange dispersing tools easily
- Digital display for precise monitoring of set and actual speeds
- Rotating knob for adjustment of speed

Ideal for the Mixing and Dispersing of Products Across a Wide Range of Viscosities  
(Up to 5000 mPas)



### High performance homogenizing.

It offers a wide speed range from 3000 – 25,000 rpm that enables users to work at high circumferential speeds even with small rotor diameters. A broad choice of dispersing elements adds versatility. Applications ranges from homogenizing waste water samples, use in laboratory reactors, dispersion tasks under vacuum / pressure and sample preparation in medical research.

### Capable of speeds between 3000 – 25,000 rpm, and volumes of 1-2000ml, this package includes everything you need to get working straight away:

- 1 x T 25 digital ULTRA-TURRAX® - Disperser
- 1 x S 25 N - 18 G - Dispersing Element Ø18 mm
- 1 x R1826 Plate stand, H800 mm
- 1 x R182 Boss head clamp, clamping range Ø6-16 mm

**Don't Miss Out on This Special Offer!**



Ref. Code:  
IKAA-T25-1024

VALID UNTIL:  
30 Nov. 2024



0800 34 24 66

sales.nz@bio-strategy.com

www.bio-strategy.com

shop.bio-strategy.com

Prices and discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices in local currency excluding GST and local delivery charges. Products advertised may be for research use only (RUO) or in vitro diagnostic use (IVD). Please refer to manufacturer's product insert for further information. For product use outside of that specified on the manufacturer's labelling and/or instructions for use, Bio-Strategy does not take responsibility.