



ESCO
WORLD CLASS. WORLDWIDE.

Airstream E Class II Microbiological Safety Cabinet

Providing operator protection, product protection and environment protection!

Only \$11,780

Usually \$16,312

Includes:
UV light and 2 power sockets

Code: ESCOAC2-4E8



External Dimensions:
1,340 x 753 x 1,400mm W x D x H
(add 57mm D for armrest) (excluding stand)

Optional adjustable
height stand available for
\$950 (usually \$1,425)
just add code ESCOSTC-4A0

- Manual sash automatically activates internal lighting and airflow when set to the marked working height
- Night set back mode reduces power consumption by a further 60%
- Multi-piece stainless steel worksurface is easy to remove for cleaning
- Ergonomic sloping front combined with tempered glass side panels allows entry of more natural light
- ULPA downflow and exhaust filters with 99.999% efficiency, 10X better than HEPA filter (99.99%) improves safety and delivers ISO Class 3 air cleanliness
- Sentinel Gold microprocessor controller with digital display shows cabinet status
- Raised armrest external to cabinet for comfortable working position
- Certified by PHE (Public Health England) for conformance to EN 12469
- Airflow sensor monitors real-time airflow for safety and alerts the user if airflow is insufficient
- Removable easy-to-clean paper catch
- Energy efficient DC ECM Motor (with 70% energy savings compared to AC motor)
- ISOCIDE™ Powder Coat – painted surfaces have silver ion impregnated coat to inhibit microbial growth
- 3-year warranty

 **bio-strategy**
Part of DKSH Group

Ref. Code:
ESCO-AC2

VALID UNTIL:
31 Dec 2024

 **DKSH**

 0800 34 24 66

 sales.nz@bio-strategy.com

 www.bio-strategy.com

 shop.bio-strategy.com

Prices and discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices in local currency excluding GST and local delivery charges. Products advertised may be for research use only (RUO) or in vitro diagnostic use (IVD). Please refer to manufacturer's product insert for further information. For product use outside of that specified on the manufacturer's labelling and/or instructions for use, Bio-Strategy does not take responsibility.