



WORLD CLASS. WORLDWIDE.

Airstream Horizontal Laminar Flow Cabinet

Why is the 1.2m Airstream Horizontal Laminar Flow Cabinet so popular?

Because it's one of the quietest and safest cabinet choices for
your sample protection requirements.

Only \$6,650

Usually \$7,671

Includes:

UV light, 2 power sockets & Front Cover

Code: ESCOLHG-4AG-F8



External Dimensions:
1,400mm (W) x 900mm (D) x 1,590mm (H)
(excluding optional stand)

Optional adjustable
height stand available for
\$950 (usually \$1,425)
just add code ESCOSTC-4AO

- ULPA filter (10x filtration efficiency than of HEPA filter) delivers ISO Class 3 air cleanliness
- Noise level: only 55.8dB(A) ; Interior lighting: 1,394 Lux
- Stainless steel spill retaining worksurface for ease of cleaning
- Sentinel™ Gold microprocessor controller with digital display showing cabinet status
- Polycarbonate fixed sash and UV resistant glass side panels allow entry of "natural" light
- Stable and self-compensating airflow
- Real time airflow monitoring with alarms
- Energy efficient DC ECM Motor (with 70% energy savings compared to AC motor)
- Standby mode provides more energy savings and maintains an ISO Class 5 work zone
- ISOCIDE™ Powder Coat – painted surfaces have silver ion impregnated coat to inhibit microbial growth
- Programmable UV timer to extend the UV lamp lifespan
- 3-year warranty



Part of DKSH Group

Ref. Code:
ESCO-LHG

VALID UNTIL:
31 Dec 2024



0800 34 24 66

sales.nz@bio-strategy.com

www.bio-strategy.com

shop.bio-strategy.com