


In Stock



Quality upright growth chamber
without humidity control

Built To Last!

Only \$ 11,975

Usually \$19,603

Code: SANYMLR-352-PE

294L with control of temperature and lighting

Includes 5 shelves and access port

- Temperature range: 0 to 50°C with lights off. +5 to 50°C with lights on
- Lighting: to 20,000 lux via 15 x 40 watt fluorescent tubes in 6 increments
- 10 program storage and each program can have up to 12-steps with repeat 1-98 times or unlimited
- Alarms for temperature (upper and lower), sensor, filter and door open
- Microprocessor PID and refrigeration capacity control with graphic LCD display
- Suitable for plant, tissue and insect culturing (insect culturing will require frequent maintenance)
- Data Logging function: up to 14 days in 6 minute intervals
- Automatic and manual defrost
- 5-year warranty



All three illuminated doors can be opened and closed for easy maintenance of the fluorescent lamps.

The structure is designed to be impenetrable to outside light.