

ESCO

WORLD CLASS. WORLDWIDE.

Ascent Opti Ductless Fume Hood



Bench mounted or
optional mobile stand

A Cost Effective Fume Hood Without Ducting!

Only \$5,995

Usually \$6,535

Code: ESCOSPD-4A1-3A1

Includes 2 x CFA-4 Code "A"
standard carbon filters

FREE (normally \$300ea)

We recommend filling out a Filtracheck form for the
Esco Technical Department to ensure your ductless
fume hood and filters meet your needs.

Contact Bio-Strategy Part of DKSH Group for the
Filtracheck form.

Offering excellent containment and ergonomic features at a great value. This compact-sized and simplified ductless fume hood is complete with internal lighting, worktop, and microprocessor control system. It also features a double-hinged sash with arm ports that provides good working flexibility and high containment levels. SPD, designed with safety and comfort in mind, also offers substantial long-term savings with its low energy consumption.

Utilising Nanocarb™ Activated Carbon Filters with optimised retention capacity, diffusion technology to ensure even filter loading and efficient perimeter clamping ensures no leakage as well as allowing for an easy filter change procedure.

- Double hinged sash mechanism with arm ports.
- Baffles for improved airflow and clearance of fumes in the work zone
- Spill retaining stainless steel work top.
- "GREEN" Environmentally friendly; does not discharge toxic gases to the environment
- Sentinel™ Silver Microprocessor Controller with digital display shows cabinet status.
- Temperature compensated airflow velocity sensor monitors airflow.
- Double hinged sash mechanism allows higher sash opening for effortless cabinet loading.
- Industry exclusive baffle design improves containment and ensures the efficient removal of chemical fumes from the work zone.
- Centrifugal fan, permanently lubricated energy efficient for reduced operating costs.
- Many optional filters available (standard, acid, mercury, sulphur, halogen, aldehyde, ammonia/amines and ethers/ chloroform.

Ref. Code: ESCO-Acent

Valid Until: 31 Mar 2025



0800 34 24 66
sales.nz@bio-strategy.dksh.com
www.bio-strategy.com
shop.bio-strategy.com



Prices and discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices in local currency excluding GST and local delivery charges. Products advertised may be for research use only (RUO) or in vitro diagnostic use (IVD). Please refer to manufacturer's product insert for further information. For product use outside of that specified on the manufacturer's labelling and/or instructions for use, Bio-Strategy does not take responsibility.