

Only \$2,725

Usually \$3,330

Code: GYROGZ-1524

Includes Snap-fit Rotor
(24 x 1.5/2.0ml)

- Max. Speed (RPM/RCF): 15,000 rpm / 21,206 xg
- Compact and light-weight design
- Quick Acc / Dec time of ≤ 15 sec
- Program memory: 100
- Selectable time mode for timer to start "AT SET SPEED" or "ON START"
- A separate "PULSE" key for a quick spin down
- "Key-lock" function to secure the parameters during operation
- 2-year warranty



FIXED ANGLE ROTOR, GRF-M2.0-24



0.2ml

Max. 16,854 xg



0.5ml

Max. 18,363 xg



1.5 / 2.0ml

Max. 21,206 xg



2.0ml with
screw cap

Max. 21,206 xg

Ref. Code: GYRO-1524

Valid Until: 28 Feb 2025