



Next Generation Moisture Analysers

With carbon fibre heating technology for fast, reliable, and accurate results

Moisture analysers are an excellent alternative to traditional oven drying methods and offer a user-friendly option compared to wet chemistry analysers. Whether you need a model for repetitive production testing or for method development in the lab, you can choose the one that best meets your needs!

Only \$2,495

Usually \$3,110

Code: OHAU30971934

Model: MB32



Perfect for routine moisture testing

- Temperature Range: 40 - 180 °C
- Drying Programs: Standard, Fast
- Shut-off Criteria: A30, A60, A90, Manual, Timed
- Method Library: 2
- Result Storage: 1
- Knob control
- 4.3" LCD display

Only \$3,650

Usually \$4,555

Code: OHAU30971935

Model: MB62



Designed for durability and flexible performance in a variety of moisture applications

- Temperature Range: 40 - 200 °C
- Drying Programs: Standard, Fast, Ramp, Step
- Shut-off Criteria: A10, A30, A60, A90, A150, Manual, Timed
- Method Library: 20
- Result Storage: 100
- Knob control
- 5" LCD display

Only \$4,150

Usually \$5,165

Code: OHAU30971936

Model: MB92



High performance moisture analysis with advanced features for the most demanding application

- Temperature Range: 40 - 200 °C
- Drying Programs: Standard, Fast, Ramp, Step
- Shut-off Criteria: A10, A30, A60, A90, A150, Manual, Timed, Free mgs
- Method Library: 20
- Result Storage: 200
- Colour easy-to-use Touchscreen display

Ref. Code: OHAU-MB

Valid Until: 31 Mar 2025



Part of DKSH Group



0800 34 24 66
sales.nz@bio-strategy.dksh.com
www.bio-strategy.com
shop.bio-strategy.com

