



Affordable 400L Laboratory Upright Refrigerator

Only \$2,245

Usually \$3,210

Code: SANYLPR-400-PE

Includes 5 height adjustable shelves

Price valid for current stock only

The LPR-400-PE Refrigerator is designed for general laboratory storage with user-friendly features.

The combination of temperature control accuracy, interior temperature uniformity, quick recovery and resistance to high ambient temperature delivers a quality storage refrigerator that characterizes PHCbi's commitment to engineering, storage safety and reliability.

- Two layer glass window in door for easy sample viewing
- Temperature Range: 4°C to 14°C selectable
- Cyclic Defrost (timer controlled)
- Low noise design; 48dB(A)
- Easy to read display with digital microprocessor control
- Forced air circulation ensures improved temperature uniformity and quick temperature recovery
- Visual & audible alarms notify end user of temperature deviation outside "safe" range
- Glass door heater prevents dew condensation without affecting the internal chamber temperature
- Access port on rear wall of cabinet for remote probe/cable access
- Outer door has cylinder keylock for security and is self-closing
- 4 castors for mobility and two levelling feet for securing in working position
- Suitable for operation on 230V 50Hz 10 amp supply
- 2-year warranty

4 Units In Stock



Ref. Code: SANY-LPR-400

Valid Until: 31 Mar 2025



0800 34 24 66
sales.nz@bio-strategy.dksh.com
www.bio-strategy.com
shop.bio-strategy.com

