



691L Biomedical Upright Freezer

Only \$9,750

Usually \$10,955

Code: SANYMDF-U731M-PE

Price valid for current stock only

**2 Units In Stock
and ready to ship**

PLUS

4 pairs of
storage bins

FREE

(valued at \$976.00)



Preservation
sub-zero



PHCbi's MDF-U731M is a large capacity laboratory freezer with high freezing power. Specifically designed to offer a larger capacity freezer with the uniformity, stability and reliability required in the life science market. Now, researchers are choosing these advanced laboratory freezers for their unrivalled performance and overall design.

- Selectable temperature range from -20°C to -30°C
- Very quiet in operation – only 42dB(A)
- Direct cooling system with manual push button defrost
- Includes 4 adjustable polyethylene coated wire shelves
- External Dimensions: 770 x 830 x 1,955mm W x D x H
- Alarms for high and low temperature, power failure, door open and compressor temperature
- Two access ports
- Castors for mobility and locating feet when in working position
- Door lock with padlock compatibility
- Dry remote alarm contacts
- Suitable for operation on 220/240V 50Hz 10 amp supply
- 5-year warranty

Ref. Code: SANY-U731M

Valid Until: 31 Mar 2025



0800 34 24 66
sales.nz@bio-strategy.dksh.com
www.bio-strategy.com
shop.bio-strategy.com



Prices and discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices in local currency excluding GST and local delivery charges. Products advertised may be for research use only (RUO) or in vitro diagnostic use (IVD). Please refer to manufacturer's product insert for further information. For product use outside of that specified on the manufacturer's labelling and/or instructions for use, Bio-Strategy does not take responsibility.