



# TwinGuard<sup>®</sup> ULT 729L Freezer

**Twice the protection  
Double your peace of mind.**

For samples you store at -86°C,  
how valuable are they?


Offering 729 litres of storage (576 x 2" cryoboxes) and enhanced safety due to the dual cooling system. Two independent refrigeration systems provide the best redundancy you can have! – if one refrigeration circuit fails – the other circuit will maintain -70°C.

- Quiet in operation – only 52dB(A)
- Filter-less condenser design: eliminates the need for routine filter cleaning
- Hot gas bypass combined with fivefold door gasket assists in minimising frosting and ice build up around door seal
- Automatic vacuum relief port for quicker access to samples after door closure
- Alarms for high and low temperature, door open and power failure
- Intelligent ECO Mode: optimising compressor operation to match the freezer's running condition and save electricity costs
- User friendly Colour Graphical Control Panel with digital display
- VIP Plus insulation: provides typically 24 % more storage capacity compared to conventionally insulated freezers with the same footprint
- Castors for mobility and locating feet when in working position
- Keyed door lock for security
- Dry remote alarm contacts
- 5-year warranty

**Only \$18,975**

**Usually \$26,500**

Code: SANYMDF-DU702VX-PE



Ref. Code: SANY-Twin729

Valid Until: 31 Mar 2025



Part of DKSH Group



0800 34 24 66  
sales.nz@bio-strategy.dksh.com  
www.bio-strategy.com  
shop.bio-strategy.com



Prices and discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices in local currency excluding GST and local delivery charges. Products advertised may be for research use only (RUO) or in vitro diagnostic use (IVD). Please refer to manufacturer's product insert for further information. For product use outside of that specified on the manufacturer's labelling and/or instructions for use, Bio-Strategy does not take responsibility.