

Smart 60L Benchtop Autoclave

Only \$35,400

Usually \$39,340

Code: PRIOB60-SMART

Included FREE

A Free Drying Kit
(valued at \$2,975)

10-year chamber warranty and 2-year warranty on balance of the autoclave

Need a flexible front loading autoclave capable of handling a variety of applications?

Complying with AS/NZS standards, the unit has a grade 316L stainless steel chamber and is inherently safe with thermal and pressure interlocks.

Suitable for waste loads, liquid/media loads, hollow glassware, tubing and porous loads.



Water Efficient



Energy Efficient



Made in Britain



Biomaster Protected



- 10 program memory
- Water cooled drain protection as standard
- Quickseal single-action door fitted with thermal and pressure locks to prevent opening at load temperature above 80°C and pressures above 0.2 bar
- Tactrol® 3 microprocessor control with simple and fully variable setting of process time and temperature with graphic indication of cycle status
- Pre-cycle vacuum in which multiple Vacuum stages are programmed with interspersed heating to achieve steam penetration
- Pulsed freesteaming to maximise efficiency of Pre-cycle vacuum phase
- Media Warming and Delayed Start functions keep autoclaved media at a usable temperature ready for immediate use
- Vacuum Assisted Cooling
- Low water level sensor
- Two loading shelves
- Forced Air Cooling
- Epoxy coated panels and frame members are treated with an anti-bacterial agent which is effective against all bacteria and fungi including MRSA
- Automatic timed free-steaming
- Pressure vessel carries a 10 year warranty and is insurance approved
- USB port, fan forced cooling and timed free-steaming
- Elements in chamber for ease of servicing & maintenance
- Many options and accessories to select from

Ref. Code: PRIO-B60

Valid Until: 31 Mar 2025

