

# DURAN® Original GL 45 Laboratory Bottles



GL 45 Laboratory Bottle, clear, with blue polypropylene screw cap and pouring ring

Packs of 10 except 5L single item

Code	Size	Was	<b>NOW</b>
SCOT218015455	1L	\$210	<b>\$189</b>
SCOT218016357	2L	\$436	<b>\$392</b>
SCOT218017353	5L	\$273	<b>\$245</b>
SCOT218012458	100ml	\$118	<b>\$106</b>
SCOT218013651	250ml	\$134	<b>\$120</b>
SCOT218014459	500ml	\$150	<b>\$135</b>

Pricing valid for current stock only

**DURAN®** Original GL laboratory bottles have become the gold standard of multifunctional laboratory bottles since their introduction in 1972.

The bottles are moulded from high-purity, clear, USP/EP Type I borosilicate glass 3.3, renowned for its excellent thermal performance, and high resistance to chemical attack.

- A combination of high-performance materials and robust ISO 4796-1 design provides a long working life
- Featuring an easy-to-read graduated volume scale, and marking field in permanent white enamel
- Permanent Retrace Code on each bottle for manufacturing lot traceability
- Supplied complete with iconic blue polypropylene cap and drip-free pouring ring for use up to 140 °C
- All components are fully autoclavable at 121 or 134 °C

Ref. Code: DWK-GL\_45

Valid Until: 30 Apr 2025