

Ascent™ Opti Formax

Formalin Dispensing
Ductless Fume Hood



Only \$8,100

Usually \$9,050

Code: ESCOSPF-3A1

Included **FREE**
1 x CFF-6 Code "F"
Aldehyde Carbon Filter

The ideal ductless fume hood for handling formalin for histological applications. The hood easily removes the irritating noxious odour and features a built-in formalin dispensing unit, carboys, a swan neck tap & foot pedal. This compact-sized and simplified ductless fume hood is complete with internal lighting, worktop, and microprocessor control system. It also features a double-hinged sash with arm ports that provides good working flexibility and high containment levels. SPF, designed with safety and comfort in mind, also offers substantial long-term savings with its low energy consumption.

Utilising a Nanocarb™ Activated Carbon Aldehyde Filter with optimised retention capacity & diffusion technology to ensure even filter loading. Efficient perimeter clamping ensures no leakage as well as allowing for an easy filter change procedure.

- Double hinged sash mechanism with arm ports
- Baffles for improved airflow and clearance of fumes in the work zone
- Low cost: no ducting required

[**Download Brochure**](#)



Comes as standard with:

- 1 x Base Cabinet with Castor Wheels
- 1 x Swan-neck Formalin Faucet
- 1 x Foot Switch
- 2 x Carboy (for Formalin)
- 1 x SS Drain Pan
- 1 x SS Worktop

Optional: Volatile Organic Compounds (VOC) Sensor
\$995 Code: VOC-SR

Provides indication of filter breakthrough warning
the user to change the filter

NB: Factory fit only – must be ordered with the hood if required

Ref. Code: ESCO-Formax

Valid Until: 31 May 2025